

# WOUND BED PREPARATION WITH PRONTOSAN®

## THE SOLUTION FOR PREVENTION AND MANAGEMENT OF BIOFILM



**1** CLEANSE AND RINSE with Prontosan® Irrigation solution



**2** WOUND DEBRIDEMENT with Prontosan® Debridement Pad



**3** REFASHION THE WOUND EDGES with Prontosan® Debridement Pad



**4** DRESS THE WOUND with Prontosan® Gel X

24/7 Action

### Steps on using the Prontosan® Debridement Pad

#### HOW TO USE FOR BEST RESULTS

#### TIP

- B. Braun branded side faces away from the wound bed
- Prontosan® Solution (as soak or rinse) can be used beforehand



**1** Open the packaging, using the integrated tray to moisten (with 15-20 mls of Prontosan Solution), covering the microfiber side of the pad.



**2** Apply light pressure, use circular or sweeping motions over areas of slough and debris. Use rounded edge for larger areas and the tapered end for difficult to access areas and shallow cavities.



**3** Following debridement, it is advisable to irrigate the wound thoroughly a second time with Prontosan Solution to cleanse the wound



**4** Apply Prontosan Wound Gel/Gel X to help prevent biofilm reformation and apply dressings as appropriate.

Cleansing Strength: One Prontosan Debridement Pad is able to remove wound debris on wounds with a surface approx, 10x 15cm (will depend on the amount of wound slough, debris etc).

### Prontosan® Debridement Pad

#### Ordering Information


















Product Size	Product Size	Article Number
Pack of 10 Pads	Droplet shape 12.76 x 9.2 cm	3908457
Pack of 3 Pads		3908456















# Prontosan® Debridement Pad

## Soft Mechanical Removal of Slough and Debris

### Prontosan® Debridement Pad field of application

WOUND DESCRIPTION / STAGE OF WOUND HEALING		LEVEL OF EXUDATE	THERAPEUTIC GOAL	PROCEDURE/ACTIONS
DEBRIDEMENT PAD	 Black necrotic wound (non-ischemic)* 	Low	Removal of dry necrotic tissue	Hydration of the wound with hydrogel. Mechanical or surgical debridement
	 Dry fibrinous wound 	Low	Removal of fibrinous tissue, Prevention of infection	Mechanical or surgical debridement if needed, wound cleansing, hydration of the wound
	 Moist fibrinous wound 	Medium to high	Removal of fibrinous tissue, Exudate management	Mechanical or surgical debridement if needed, wound cleansing, enhancing moist environment to promote autolytic debridement
	  Fibrinous to granulating wound 	Medium to high	Removal of fibrinous tissue, Exudate management, Promotion of granulation	Wound cleansing, enhancing moist environment to promote autolytic debridement
	 Granulating wound 	Medium	Exudate management, Promotion of granulation	Wound cleansing, moist environment, protection of newly formed tissue and of surrounding skin
	 Granulating to epithelialising wound 	Low	Absorption of remaining exudate, Promotion of epithelialisation	Wound cleansing, protection of newly formed tissue and of surrounding skin
	 Epithelialising wound 	None	Promotion of epithelialisation	Protection of newly formed tissue
	 Infected wound 	Medium to high	Reduction of bacterial load, Exudate management, Odour control	Wound cleansing, prevent contamination from external bacteria, enhance moist wound environment

### Application protocol with Prontosan® Debridement Pad

CHRONIC WOUND - GRANULATING		Soak with Prontosan Irrigation Solution	Debride with Pad & Solution	Consider Prontosan Wound Gel X
<ul style="list-style-type: none"> <li>High risk patient*</li> <li>Low exudate</li> </ul> 	<ul style="list-style-type: none"> <li>Cleans</li> <li>Prevents biofilm/complications</li> </ul>			
<b>CHRONIC WOUND</b>				
<ul style="list-style-type: none"> <li>Light slough</li> <li>Low exudate</li> </ul> 	<ul style="list-style-type: none"> <li>Cleans</li> <li>Prevents biofilm/complications</li> </ul>			
<b>CHRONIC WOUND - CRITICALLY COLONISED/INFECTED</b>				
<ul style="list-style-type: none"> <li>Medium/high exudate</li> <li>Static wound</li> <li>Slough</li> </ul> 	<ul style="list-style-type: none"> <li>Cleans</li> <li>Prevents biofilm/complications</li> </ul>			

\*High risk patient: Co-morbidities such as Diabetes, immuno-compromised, steroidal use, patients with previous wound infection and or biofilm and slough.