

Metal handle – a better grip than metal?

Comparison SQ.line® mallet vs. Aesculap® silicone version



Aesculap® SQ.line® Aesculap® Silicone Version



Key Facts at a glance:

Comparing grip forces, the **SQ.line® mallet** offers a **safer grip** than the silicone version as follows:

- **Safer grip of SQ.line®...**
...regarding distal end (pull force): **>10%**
...regarding proximal end (push force): **70%**.
- Against **rotation** (torsion force) the metal and silicon version are **at the same level**.
- There are improvements concerning grip properties of the metal version.

- The innovative design of the SQ.line® handles - new handle shape, contours and textural features - requires less gripping force to achieve good control during use. This translates to **less fatigue** on the operator!



Executive Summary - Test performance:


Characteristics of testing persons:

10 testing persons (5 male , 5 female, age 20-54, glove size S+M+L)

Test 1: Push & Pull Forces

The different mallets, with handle designs in "metal" and "silicone", were fixed in a machine which measures axial push and pull forces. For this test, each handle was lubed with Vaseline to simulate the slippery condition of surgical use.

The force in Newton (N) was tested and recorded for the point at which the mallet began to slip within the hand of the testing person.



Test results:

	Ø Pull force	Ø Push force
Aesculap® silicone handle	74 N (100%) (Deviation: 22N-108N)	55 N (100%) (Deviation: 30N-78N)
SQ.line® metal handle	83 N (113%) (Deviation: 56N-127N)	94 N (170%) (Deviation: 66N-152N)

Test 2: Torsion force

The different mallets were fixed in a machine which measures torsion forces.

Again the handles were lubed.

For each handle, the force in newton meter (nm) was tested and recorded for the point at which the mallet starts to slip in the hand of the testing person.



Test results:

	Ø Torsion force
Aesculap® silicone handle	1.702 Nm (100%) (Deviation: 1.30nm-2.24nm)
SQ.line® metal handle	1.756 Nm (103%) (Deviation: 0.97nm-2,47nm)

AESCULAP® – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany

Phone +49 7461 95-0 | Fax +49 7461 95-26 00 | www.aesculap.com

The main product trademark 'Aesculap' and the product trademark 'SQ.line' are registered trademarks of Aesculap AG. Subject to technical changes. All rights reserved. This datasheet may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

D-XX19046 0520/PDF/1