

PRODUCTS WITH
NRFit™
CONNECTORS



Reducing Risks of Misconnections in Regional Anaesthesia¹

Connecting what belongs together

At B. Braun, we are ready for the change.

As your industry partner we're here to help you stay on top of evolving industry safety standards and transition to products with NRFit® connectors.

Helping equip health care facilities with products with NRFit® connectors is part of B. Braun's commitment to improve patient and facility outcomes.

Visit: <https://www.bbraun.com.au/en/products-and-therapies/pain-therapy/nrfit> to find out more.



To improve patient safety by reducing the risk of inadvertent misconnections, Luer medical device connectors for neuraxial applications are changing to meet international design standards – NRFit®.²

NRFit® connectors are 20% smaller in diameter than Luer connectors, preventing medical devices meant for neuraxial administration from connecting to devices used for IV, enteral and other applications and vice versa.²

NRFit® Connector



Luer Connector

Luer



Epidural Needle



Spinal Needle



LOR Syringe



Procedural Syringe (Lock)



Procedural Syringe (Slip)



Catheter Connector



Epidural Flat Filter

NRFit™



B. Braun Products with NRFit™ Connectors

Visible Improvements

- NRFit® connector is longer than Luer connector
- 20% smaller connector diameter
- Yellow product colouring or yellow cap where feasible
- Collared Slip syringes
- Addition of prism in spinal needle hub for CSF identification
- NRFit® needle stylets are a transparent colour rather than opaque as seen in Luer needles
- NRFit® product codes include letter "N" where possible
- Product Labeling denotes NRFit®

Excluded Products

- Spinal introducer needle
- Hypodermic needles

References

1. Dedicated connectors intended to help reduce the risk of misconnections. Yellow colour coding where feasible allows faster identification of Regional Anaesthesia devices.

2. ISO 80369-6:2016 - small bore connectors for liquids and gases in healthcare applications - Part 6: connectors for neuraxial applications.

B. Braun Australia Pty Ltd | Level 5, 7-9 Irvine Place, Bella Vista NSW 2153,
Australia | Tel: 1800 251 705 | info.au@bbraun.com | www.bbraun.com.au

B. Braun New Zealand | PO Box 37353 Parnell, Auckland 1052, New Zealand |
Customer Care 0800 227 286 | Fax (09) 373 5601 | www.bbraun.co.nz

BANZ HC E0778 02/22 | © 2022 B. Braun Australia