

REUSABLE HANDLE

for use with monopolar single use shafts



	Description	Art. No.	Pcs. per pack
 1:1	Handle without ratchet	P0958R	1
 1:1	Handle with ratchet, ratchet can be deactivated	P0959R	1

AESCULAP® – a B. Braun brand

B. Braun Australia Pty Ltd | Level 2, 4 Burbank Place, Norwest NSW 2153, Australia | Tel.1800 251 705 | info.au@bbraun.com | www.bbraun.com.au

©2025 B. Braun Australia Pty Ltd. ABN 56 002 945 155. All rights reserved. Aesculap is a registered trademark of B. Braun, Australia. This brochure is intended for use by healthcare professionals only. It is provided for informational purposes and does not constitute medical advice. Always refer to the product's Instructions for Use and consult a qualified healthcare provider for clinical decisions.

AESCULAP®

AdTec® single use

Combination of sterile single use shafts and reusable monopolar handles



AESCULAP® AdTec® single use



MONOPOLAR SINGLE USE SHAFTS

- Combination of sterile single use shafts and reusable monopolar handles
- „One-Click“ design:
 - assembly and disassembly in one step only
 - easier handling
- Insulated up to the jaws (double action jaws)
- Easier force adjustment through mechanical force transmission from the jaw part to the handle
- Shaft 360° rotatable
- Economic price performance ratio
- Pack containing 10 pieces (single packed)

SCISSORS

Reliability: good cutting properties.¹



	Description	Art. No.	Length	Diameter	Pcs. per pack
	Scissors, tips slight curved with fine tip	PO885SU	310 mm	5 mm	10
	Hook scissors	PO886	310 mm	5 mm	10
	Mini-METZENBAUM scissors, tips curved to the left	PO887	310 mm	5 mm	10
	METZENBAUM scissors, tips curved to the left	PO888	310 mm	5 mm	10
	METZENBAUM scissors, tips curved to the left	PO889	420 mm	5 mm	10

FORCEPS

Reliability: precise grasping and dissection of tissue.¹



	Description	Art. No.	Length	Diameter	Pcs. per pack
	MARYLAND fixation and dissecting forceps	PO891SU	310 mm	5 mm	10
	Pike-mouth forceps	PO892SU	310 mm	5 mm	10
	Atraumatic universal grasping forceps, fenestrated	PO893SU	310 mm	5 mm	10
	BABCOCK tissue – grasping forceps	PO894SU	310 mm	5 mm	10

¹ According to feedback within the scope of a survey of 76 surgeons that use monopolar or bipolar AdTec instruments. The respective survey was conducted by Aesculap.