



AESCULAP® Trocar Solution

Single Use System

AESCULAP® Single Use Trocar Solution

Optical – Dilating – Bladed / Shielded – HASSON – Sleeves



Clear tip allows visualization of abdominal wall layers

- Optical balloon trocars with stopcock and ring / silicone cone
- Optical trocars with ribbed surface and stopcock
- 11 mm / 12 mm trocars with universal seal and removable trocar head
- Firm hold in the abdominal wall*




Optical Balloon Trocars

Optical Trocars

Optical Trocar Sets

Article Number	Description	Diameter (mm)	Length (mm)	Packing unit
Optical Balloon Trocars				
EK500SU	Optical balloon trocar (ring)	5	70	1 Pack of 5 pieces
EK501SU	Optical balloon trocar (ring)	5	100	1 Pack of 5 pieces
EK502SU	Optical balloon trocar with universal seal (ring)	11	100	1 Pack of 5 pieces
EK503SU	Optical balloon trocar with universal seal (ring)	12	100	1 Pack of 5 pieces
EK504SU	Optical balloon trocar with universal seal (silicone cone)	12	100	1 Pack of 5 pieces



Optical Trocars

EK530SU	Optical trocar	5	70	1 Pack of 5 pieces
EK531SU	Optical trocar	5	100	1 Pack of 5 pieces
EK534SU	Optical trocar with universal seal	11	100	1 Pack of 5 pieces
EK536SU	Optical trocar with universal seal	12	75	1 Pack of 5 pieces
EK537SU	Optical trocar with universal seal	12	100	1 Pack of 5 pieces
EK538SU	Optical trocar with universal seal	12	150	1 Pack of 5 pieces

Optical Trocar Sets

(including 2 trocar sleeves and 1 trocar pin)

EK540SU	Optical trocar set	5	70	1 Pack of 5 pieces
EK541SU	Optical trocar set	5	100	1 Pack of 5 pieces
EK542SU	Optical trocar set with universal seal	11	100	1 Pack of 5 pieces
EK543SU	Optical trocar set with universal seal	12	100	1 Pack of 5 pieces

AESCULAP® Single Use Trocar Solution

Optical – Dilating – Bladed / Shielded – HASSON – Sleeves



Dilating trocars with lateral wings

- Dilating balloon trocars with stopcock and ring / silicone cone
- Dilating trocars with ribbed surface and stopcock
- 11 mm / 12 mm / 15 mm trocars with universal seal and removable trocar head
- Firm hold in the abdominal wall*




Dilating Balloon Trocars

Dilating Trocars

Dilating Trocar Sets

Article Number	Description	Diameter (mm)	Length (mm)	Packing unit
Dilating Balloon Trocars				
EK510SU	Dilating balloon trocar (ring)	5	70	1 Pack of 5 pieces
EK511SU	Dilating balloon trocar (ring)	5	100	1 Pack of 5 pieces
EK512SU	Dilating balloon trocar with universal seal (ring)	11	100	1 Pack of 5 pieces
EK513SU	Dilating balloon trocar with universal seal (ring)	12	100	1 Pack of 5 pieces
EK514SU	Dilating balloon trocar with universal seal (silicone cone)	12	100	1 Pack of 5 pieces



Dilating Trocars

EK570SU	Dilating trocar	5	70	1 Pack of 5 pieces
EK571SU	Dilating trocar	5	100	1 Pack of 5 pieces
EK574SU	Dilating trocar with universal seal	11	100	1 Pack of 5 pieces
EK577SU	Dilating trocar with universal seal	12	100	1 Pack of 5 pieces
EK578SU	Dilating trocar with universal seal	12	150	1 Pack of 5 pieces
EK579SU	Dilating trocar with universal seal	15	100	1 Pack of 5 pieces

Dilating Trocar Sets

(including 2 trocar sleeves and 1 trocar pin)

EK580SU	Dilating trocar set	5	70	1 Pack of 5 pieces
EK581SU	Dilating trocar set	5	100	1 Pack of 5 pieces
EK582SU	Dilating trocar set with universal seal	11	100	1 Pack of 5 pieces
EK583SU	Dilating trocar set with universal seal	12	100	1 Pack of 5 pieces

AESCULAP® Single Use Trocar Solution

Optical – Dilating – Bladed / Shielded – HASSON – Sleeves



Blue activation button for safety shield

- Bladed / shielded balloon trocars with stopcock and ring / silicone cone
- Bladed / shielded trocars with ribbed surface and stopcock
- 11 mm / 12 mm / 15 mm trocars with universal seal and removable trocar head
- Firm hold in the abdominal wall*



Bladed / Shielded
Balloon Trocars

Bladed / Shielded
Trocars

Bladed / Shielded
Trocar Sets

Article Number	Description	Diameter (mm)	Length (mm)	Packing unit
Bladed / Shielded Balloon Trocars				
EK505SU	Bladed / shielded balloon trocar (ring)	5	70	1 Pack of 5 pieces
EK506SU	Bladed / shielded balloon trocar (ring)	5	100	1 Pack of 5 pieces
EK507SU	Bladed / shielded balloon trocar with universal seal (ring)	11	100	1 Pack of 5 pieces
EK508SU	Bladed / shielded balloon trocar with universal seal (ring)	12	100	1 Pack of 5 pieces
EK509SU	Bladed / shielded balloon trocar with universal seal (silicone cone)	12	100	1 Pack of 5 pieces



Silicone cone



Ring

Bladed / Shielded Trocars

EK550SU	Bladed / shielded trocar	5	70	1 Pack of 5 pieces
EK551SU	Bladed / shielded trocar	5	100	1 Pack of 5 pieces
EK554SU	Bladed / shielded trocar with universal seal	11	100	1 Pack of 5 pieces
EK557SU	Bladed / shielded trocar with universal seal	12	100	1 Pack of 5 pieces
EK559SU	Bladed / shielded trocar with universal seal	15	100	1 Pack of 5 pieces

Bladed / Shielded Trocar Sets

(including 2 trocar sleeves and 1 trocar pin)

EK560SU	Bladed / shielded trocar set	5	70	1 Pack of 5 pieces
EK561SU	Bladed / shielded trocar set	5	100	1 Pack of 5 pieces
EK562SU	Bladed / shielded trocar with universal seal	11	100	1 Pack of 5 pieces
EK563SU	Bladed / shielded trocar with universal seal	12	100	1 Pack of 5 pieces

AESCULAP® Single Use Trocar Solution

Optical – Dilating – Bladed / Shielded – HASSON – Sleeves



HASSON Trocars in different options

- HASSON balloon trocars with stopcock and ring / silicone cone
- HASSON trocar with stopcock and cone
- 11 mm / 12 mm trocars with universal seal and removable trocar head



HASSON Balloon Trocars



HASSON Trocar

Article Number	Description	Diameter (mm)	Length (mm)	Packing unit
HASSON Balloon Trocars				
EK515SU	HASSON balloon trocar with universal seal (ring)	11	100	1 Pack of 5 pieces
EK516SU	HASSON balloon trocar with universal seal (ring)	12	100	1 Pack of 5 pieces
EK517SU	HASSON balloon trocar with universal seal (silicone cone)	12	100	1 Pack of 5 pieces
HASSON Trocar				
EK588SU	HASSON trocar with universal seal and cone	11	100	1 Pack of 5 pieces
EK589SU	HASSON trocar with universal seal and cone	12	100	1 Pack of 5 pieces



AESCULAP® Single Use Trocar Solution

Optical – Dilating – Bladed / Shielded – HASSON – Sleeves



Various trocar sleeves single packed


- Balloon trocar sleeves with stopcock and ring / silicone cone
- Trocar sleeves with ribbed surface and stopcock
- 11 mm / 12 mm trocar sleeves with universal seal and removable trocar head
- Firm hold in the abdominal wall*



Balloon Trocar Sleeves

Trocar Sleeves

Article Number	Description	Diameter (mm)	Length (mm)	Packing unit
Balloon Trocar Sleeves				
EK520SU	Trocar sleeve with balloon (ring)	5	70	1 Pack of 5 pieces
EK521SU	Trocar sleeve with balloon (ring)	5	100	1 Pack of 5 pieces
EK522SU	Trocar sleeve with balloon and universal seal (ring)	11	100	1 Pack of 5 pieces
EK523SU	Trocar sleeve with balloon and universal seal (ring)	12	100	1 Pack of 5 pieces
EK524SU	Trocar sleeve with balloon and universal seal (silicone cone)	12	100	1 Pack of 5 pieces



Trocar Sleeves				
EK590SU	Trocar sleeve	5	70	1 Pack of 5 pieces
EK591SU	Trocar sleeve	5	100	1 Pack of 5 pieces
EK594SU	Trocar sleeve with universal seal	11	100	1 Pack of 5 pieces
EK597SU	Trocar sleeve with universal seal	12	100	1 Pack of 5 pieces
EK598SU	Trocar sleeve with universal seal	12	150	1 Pack of 5 pieces

*based on customer feedback in the context of the usability evaluation. The fixation in the abdominal wall of animal specimen was rated firm by 9 of 11 customers.

AESCULAP® – a B. Braun brand

B. Braun Australia Pty Ltd | Norwest Business Park NSW, 2153 | Tel. 1800 251 705 | info.au@bbraun.com | bbraun.com.au

B. Braun New Zealand Pty Ltd | Level 8, 139 Quay Street, Auckland, 1010, New Zealand | Tel. 0800 227 286 | info.nz@bbraun.com | bbraun.co.nz

Manufacturer acc. to MDD 93/42: Aesculap AG, Am Aesculap-Platz, 78532 Tuttlingen, Germany