



AdTec[®] Monopolar

Laparoscopic, Monopolar,
Reusable Instruments


AdTec® Monopolar



For years, AESCULAP® has been actively involved in the development of laparoscopic techniques in the instruments sector, introducing products such as AdTec® Single use, AdTec® needle holders as well as the Multi-Fire Clip Applier Challenger Ti-P. Here, great value has always been placed on quality, user friendliness, efficiency and customer satisfaction.

In addition, AESCULAP® offers a comprehensive package of services ranging from planning and consultancy, instrument screening as well as training and qualification up to Instrument repair.

Why not enjoy the benefits of the reusable monopolar instrument line?



One-Click Design

- Assembly and disassembly in one step only
- Easier handling

Ratchet Mechanism


- Ratchet can be deactivated
- Easier handling



Jaw Mechanism

- Non-protruding joints during opening and closure of jaw Parts (for the majority of instruments, for details see item list)
- Limits tissue or suture material that can get caught in the jaw mechanism (for the majority of instruments, for details see item list)

AdTec® Monopolar



Peek Insulation

- Enhanced spark protection and resisting characteristic
- Quality instruments specifically designed for CSSD processing

Economical Price Performance Ratio

- Enhanced modular design aids financial benefits

Universal


- Modular design consisting of four Components, facilitating the disassembly, reprocessing, maintenance and sterile supply











Sterilisation

- Sterilisation of the assembled Instrument possible

AdTec® Monopolar


Scissors










	Complete Instrument		Single Parts			Individual Components	
			Jaw insert	Insulated outer tube	Handle		
 1/1	METZENBAUM Scissors, serrated, tips curved to the left, with carbide inserts	5 mm / 220 mm	P0631R 	PM970R 	P0958R 	PM971R  Ø 5/10 mm	P0958R 
 1/1	METZENBAUM Scissors, serrated, tips curved to the left, with carbide inserts	5 mm / 310 mm	P0603R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm	P0959R 
 1/1	METZENBAUM Scissors, serrated, tips curved to the left, with carbide inserts	5 mm / 420 mm	P0746R 	PM981R 	P0958R 		
 1/1	METZENBAUM Scissors, serrated, tips curved to left, with carbide inserts, PEEK insulation of the distal shaft	5 mm / 310 mm	PM697R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm	PM953R 
 1/1	CROCE/OLMI Dissecting and suture scissors with knot pusher	5 mm / 310 mm	P0727R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm	PM954R 
 1/1	METZENBAUM Scissors, serrated, tips curved to the left, with carbide inserts	10 mm / 310 mm	PM604R 	PM974R 	P0958R 		
 1/1	Noir® METZENBAUM Scissors, serrated, tips curved to the left, with carbide inserts	5 mm / 310 mm	P0790R 	PM973R 	P0958R 		P0957R 
 1/1	Noir® Mini-METZENBAUM Scissors, serrated, tips curved to the left, with carbide inserts	5 mm / 310 mm	P0791R 	PM973R 	P0958R 		

AdTec® Monopolar


Scissors



	Complete Instrument		Single Parts			Individual Components	
			Jaw insert	Insulated outer tube	Handle		
 18 mm 1/1	PO030R Medium-METZENBAUM Scissors, serrated, tips curved to the left, with carbide inserts	5 mm / 310 mm	P0715R 	PM970R 	P0958R 		P0958R 
 18 mm 1/1	PO031R Medium-METZENBAUM Scissors, serrated, tips curved to the left, with carbide inserts*	5 mm / 420 mm	P0716R 	PM981R 	P0958R 		P0959R 
 14 mm 1/1	PO023R Mini-METZENBAUM Scissors, serrated, tips curved to the left, with carbide inserts	5 mm / 220 mm	P0738R 	PM970R 	P0958R 	PM971R  Ø 5/10 mm	P0959R 
 14 mm 1/1	PO011R Mini-METZENBAUM Scissors, serrated, tips curved to the left, with carbide inserts	5 mm / 310 mm	P0655R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm	PM953R 
 14 mm 1/1	PO028R Mini-METZENBAUM Scissors, serrated, tips curved to the left, with carbide inserts	5 mm / 420 mm	P0754R 	PM981R 	P0958R 		PM954R 
 13 mm 1/1	Mini-METZENBAUM Scissors, serrated, tips curved to the left, with carbide inserts, PEEK insulation of the distal shaft	10 mm / 310 mm	PM698R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm	PM954R 
 13 mm 1/1	Mini-METZENBAUM Scissors, serrated, tips curved to the left, with carbide inserts, PEEK insulation of the distal shaft	5 mm / 310 mm	PM710R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm	PM957R 

AdTec® Monopolar


Scissors











		Complete Instrument	Single Parts			Individual Components	
			Jaw insert	Insulated outer tube	Handle		
	PO001R	Hook scissors, blunt	5 mm / 310 mm	P0600R	PM973R	P0958R	PM975R
	1/1						Ø 5/10 mm
	PO025R	Hook scissors, blunt	5 mm / 420 mm	P0744R	PM981R	P0958R	
	1/1						
	PO002R	Peritoneal scissors, serrated	5 mm / 310 mm	P0601R	PM973R	P0958R	PM975R
	1/1						Ø 5/10 mm
	PO006R	Peritoneal scissors with hook, serrated	5 mm / 310 mm	P0625R	PM973R	P0958R	PM975R
	1/1						Ø 5/10 mm
		Micro-dissecting scissors, tips curved to the left	5 mm / 220 mm	PM630R	PM970R	P0958R	PM971R
	1/1						Ø 5/10 mm
	PO003R	Micro-dissecting scissors, tips curved to the left	5 mm / 310 mm	P0602R	PM973R	P0958R	PM975R
	1/1						Ø 5/10 mm
	PO014R	POTTS de MARTELL Mod. SAID DUROTIP Endo vessel scissors	10 mm / 310 mm	P0682R	PM992R	P0958R	
	1/1	Instrument with protruding jaw Parts					

AdTec® Monopolar


Dissecting Forceps













	Complete Instrument		Single Parts			Individual Components	
			Jaw insert	Insulated outer tube	Handle		
 1/1	MARYLAND Fixation and dissecting forceps, curved	5 mm / 220 mm	P0632R 	PM970R 	P0958R 	PM971R  Ø 5/10 mm	P0958R 
 1/1	MARYLAND Fixation and dissecting forceps, curved	5 mm / 310 mm	P0608R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm	P0959R 
 1/1	MARYLAND Fixation and dissecting forceps, curved	5 mm / 420 mm	P0752R 	PM981R 	P0958R 		
 1/1	MARYLAND Fixation and dissecting forceps, insulation of the distal shaft Instrument with protruding jaw Parts	5 mm / 310 mm	PM699R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm	PM953R 
 1/1	MARYLAND Fixation and dissecting forceps, curved with cross-teething	5 mm / 310 mm	P0717R 	PM973R 	P0958R 		PM954R 
 1/1	MARYLAND Fixation and dissecting forceps, curved with cross-teething	5 mm / 420 mm	P0718R 	PM981R 	P0958R 		PM957R 

AdTec® Monopolar


Dissecting Forceps












		Complete Instrument	Single Parts			Individual Components		
			Jaw insert	Insulated outer tube	Handle			
 20 mm 1/1	PO195R	CROCE/OLMI Fixation and dissection forceps, curved, jaw Parts fenestrated	5 mm / 310 mm	 PO740R	 PM973R	 P0958R	 PM975R Ø 5/10 mm	 P0958R
 26 mm 1/1	PO100R	Forceps with atraumatic teeth	5 mm / 310 mm	 PO605R	 PM973R	 P0958R	 PM975R Ø 5/10 mm	 P0959R
 19 mm 1/1		Fixation and dissecting forceps, slim	5 mm / 310 mm	 PM614R	 PM973R	 P0958R	 PM975R Ø 5/10 mm	 P0958R
 25 mm 1/1		Instrument with protruding jaw Parts		 PM647R	 PM973R	 P0958R	 PM975R Ø 5/10 mm	 PM953R
 16 mm 1/1	PO187R	PSALMON Precise-dissector, 90°	5 mm / 310 mm	 PO730R	 PM973R	 P0958R	 PM975R Ø 5/10 mm	 PM954R
 27 mm 1/1	PO231R	MIXTER Dissecting clamp, 90°	5 mm / 310 mm	 PO748R	 PM973R	 P0958R	 PM975R Ø 5/10 mm	 PM957R

AdTec® Monopolar


Dissecting Forceps













	Complete Instrument		Single Parts			Individual Components	
			Jaw insert	Insulated outer tube	Handle		
 28 mm 1/1	PO103R	OVERHOLT Dissecting clamp, 90° 10 mm / 310 mm	PO606R 	PM975R 	P0958R 	PM973R  Ø 5 mm	P0958R 
 27 mm 1/1		OVERHOLT Dissecting clamp, 90°, PEEK insulation of the distal shaft 10 mm / 370 mm Instrument with protruding jaw Parts	PM701R 	PM978R 	P0958R 	PM976R  Ø 5 mm	P0959R 
 27 mm 1/1	PO219R	OVERHOLT Dissecting forceps with fine tip, 90° 10 mm / 370 mm	P0772R 	PM978R 	P0958R 	PM976R  Ø 5 mm	
 40 mm 1/1	PO141R	OVERHOLT Dissecting clamp, 60° 10 mm / 310 mm	PO661R 	PM975R 	P0958R 	PM973R  Ø 5 mm	PM953R 
 39 mm 1/1		OVERHOLT Dissecting clamp, 60°, PEEK insulation of the distal shaft 10 mm / 310 mm Instrument with protruding jaw Parts	PM700R 	PM975R 	P0958R 	PM973R  Ø 5 mm	PM954R 
 40 mm 1/1	PO220R	OVERHOLT Dissecting forceps with fine tip, 60° 10 mm / 370 mm	P0773R 	PM978R 	P0958R 	PM976R  Ø 5 mm	
 30 mm 1/1	PO218R	CRILE Fixation and dissecting forceps, curved 5 mm / 310 mm Instrument with protruding jaw Parts	P0771R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm	PM957R 
 30 mm 1/1	PO213R	CRILE Fixation and dissecting forceps, curved 5 mm / 370 mm Instrument with protruding jaw Parts	P0766R 	PM976R 	P0958R 	PM978R  Ø 5/10 mm	

AdTec® Monopolar


Dissecting Forceps








	Complete Instrument	Single Parts			Individual Components
		Jaw insert	Insulated outer tube	Handle	
 17 mm 1/1	PO147R Fixation forceps with jaw groove, slim tapering jaw tip 5 mm / 310 mm	PO678R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm 
 18 mm 1/1	Atraumatic fixation and dissecting forceps 5 mm / 220 mm	PO633R 	PM970R 	P0958R 	PM971R  Ø 5/10 mm 
 18 mm 1/1	Atraumatic fixation and dissecting forceps 5 mm / 310 mm	PO607R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm 
 18 mm 1/1	Atraumatic fixation and dissecting forceps 5 mm / 420 mm	PO747R 	PM981R 	P0958R 	
 19 mm 1/1	Atraumatic universal forceps 5 mm / 370 mm	PO767R 	PM976R 	P0958R 	PM978R  Ø 5/10 mm 
 16 mm 1/1	Delicate dissecting forceps 5 mm / 310 mm	PO667R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm 
 17 mm 1/1	Fixation and dissecting forceps with jaw groove 5 mm / 310 mm	PO735R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm
 17 mm 1/1	Fixation and dissecting forceps with jaw groove 5 mm / 420 mm	PO753R 	PM981R 	P0958R 	PM975R  Ø 5/10 mm

AdTec® Monopolar

Grasping Forceps










	Complete Instrument	Single Parts			Individual Components
		Jaw insert	Insulated outer tube	Handle	
 20 mm 1/1	P0109R Crocodile forceps, extra-long jaws 5 mm / 310 mm	 P0615R	 PM973R	 P0959R	 PM975R Ø 5/10 mm  P0958R
 22 mm 1/1	P0127R Universal fixation forceps with jaw grooves 5 mm / 310 mm	 P0644R	 PM973R	 P0959R	 PM975R Ø 5/10 mm  P0959R
 20 mm 1/1	P0221R Universal grasping forceps with reservoir in the jaws 5 mm / 310 mm	 P0774R	 PM973R	 P0959R	 PM975R Ø 5/10 mm  PM953R
 20 mm 1/1	P0230R Universal grasping forceps with reservoir in the jaws 5 mm / 420 mm	 P0783R	 PM981R	 P0959R	 PM953R  PM954R  PM957R

Grasping Forceps

AdTec® Monopolar

Grasping Forceps




	Complete Instrument		Single Parts			Individual Components		
			Jaw insert	Insulated outer tube	Handle			
 28 mm 1/1	PO233R	Atraumatic universal forceps with grooves	5 mm / 310 mm	P0786R	PM973R	P0959R	PM975R Ø 5/10 mm	P0958R
 28 mm 1/1	PO234R	Atraumatic universal forceps with grooves	5 mm / 420 mm	P0787R	PM981R	P0959R	PM975R Ø 5/10 mm	P0959R
 27 mm 1/1	PO235R	Atraumatic universal forceps with grooves	5 mm / 310 mm	P0788R	PM973R	P0959R	PM975R Ø 5/10 mm	
 27 mm 1/1	PO236R	Atraumatic universal forceps with grooves	5 mm / 420 mm	P0789R	PM981R	P0959R		PM953R
 25 mm 1/1	PO216R	Atraumatic grasping forceps, jaw part fenestrated	5 mm / 310 mm	P0769R	PM973R	P0958R	PM975R Ø 5/10 mm	PM954R
 25 mm 1/1	PO217R	Atraumatic grasping forceps, jaw part fenestrated	5 mm / 420 mm	P0770R	PM981R	P0958R		
 24 mm 1/1	PO128R	Atraumatic universal grasping forceps, jaw part fenestrated	5 mm / 310 mm	P0647R	PM973R	P0959R	PM975R Ø 5/10 mm	PM957R
 24 mm 1/1	PO229R	Atraumatic universal grasping forceps, jaw part fenestrated	5 mm / 420 mm	P0782R	PM981R	P0959R	PM975R Ø 5/10 mm	

Grasping Forceps

AdTec® Monopolar


Grasping Forceps





















	Complete Instrument		Single Parts			Individual Components		
			Jaw insert	Insulated outer tube	Handle			
 30 mm 1/1	PO191R	Universal fixation and dissecting forceps, jaw part fenestrated	5 mm / 310 mm	P0736R	PM973R	P0959R	PM975R Ø 5/10 mm	P0958R
 30 mm 1/1	PO185R	CROCE/OLMI Atraumatic grasping forceps	5 mm / 310 mm	P0728R	PM973R	P0959R	PM975R Ø 5/10 mm	P0959R
 30 mm 1/1	PO232R	CROCE/OLMI Atraumatic grasping forceps, jaw part fenestrated	5 mm / 310 mm	P0785R	PM973R	P0959R	PM975R Ø 5/10 mm	P0959R
 19 mm 1/1	PO222R	Short universal fixation and dissecting forceps, jaw part fenestrated	5 mm / 310 mm	P0775R	PM973R	P0959R	PM975R Ø 5/10 mm	PM953R
 19 mm 1/1	PO223R	Short universal fixation and dissecting forceps, jaw part fenestrated	5 mm / 420 mm	P0776R	PM981R	P0959R		PM954R
 35 mm 1/1	PO112R	DORSEY Intestinal grasping forceps	5 mm / 310 mm	P0618R	PM973R	P0959R	PM975R Ø 5/10 mm	PM957R
 35 mm 1/1		DORSEY Intestinal grasping forceps	5 mm / 370 mm	PM618R	PM976R	P0959R	PM978R Ø 5/10 mm	
 35 mm 1/1	PO204R	DORSEY Intestinal grasping forceps	5 mm / 420 mm	P0755R	PM981R	P0959R		
		Instrument with protruding jaw Parts						

AdTec® Monopolar

Grasping Forceps

















	Complete Instrument		Single Parts			Individual Components	
			Jaw insert	Insulated outer tube	Handle		
 PO113R 1/1	BABCOCK Tissue-grasping forceps	5 mm / 310 mm	PO619R 	PM973R 	PO959R 	PM975R 	PO958R 
 1/1	BABCOCK Tissue-grasping forceps For use with flexible 13 mm trocars only	5 mm / 310 mm	PM640R 	PM973R 	PO959R 	PM975R 	PO959R 
 PO114R 1/1	BABCOCK Tissue-grasping forceps	10 mm / 310 mm	PO620R 	PM974R 	PO959R 		
 PO205R 1/1	BABCOCK Tissue-grasping forceps	10 mm / 420 mm	PO757R 	PM982R 	PO959R 		PM953R 
 PO212R 1/1	BABCOCK Tissue-grasping forceps	5 mm / 310 mm	PO765R 	PM973R 	PO959R 	PM975R 	PM954R 
 PO209R 1/1	BABCOCK Tissue-grasping forceps	5 mm / 420 mm	PO762R 	PM981R 	PO959R 		
 PO196R 1/1	BABCOCK Tissue-grasping forceps	10 mm / 370 mm	PO741R 	PM977R 	PO959R 		
 PO197R 3/4	MEYER Tissue grasping and manipulation forceps, jaw length 90 mm	10 mm / 370 mm	PO742R 	PM977R 	PO959R 		PM957R 

Grasping Forceps

AdTec® Monopolar

Grasping Forceps


















	Complete Instrument		Single Parts			Individual Components	
			Jaw insert	Insulated outer tube	Handle		
 35 mm 1/1	PO104R ALLIS Atraumatic tissue grasping forceps	5 mm / 310 mm	PO610R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm	P0958R 
 30 mm 1/1	ALLIS Atraumatic grasping forceps	5 mm / 310 mm	PM670R 	PM973R 	P0959R 	PM975R  Ø 5/10 mm	P0959R 
 36 mm 1/1	LINDER Endo lung grasping forceps	5 mm / 220 mm	PM696R 	PM970R 	P0958R 	PM971R  Ø 5/10 mm	
 36 mm 1/1	LINDER Endo lung grasping forceps	5 mm / 310 mm	PM687R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm	PM953R 
 29 mm 1/1	DUVAL Atraumatic endo parenchyma fixation forceps	5 mm / 310 mm	PM609R 	PM973R 	P0959R 	PM975R  Ø 5/10 mm	PM954R 
 32 mm 1/1	DUVAL Atraumatic endo parenchyma fixation forceps	10 mm / 370 mm	PM612R 	PM977R 	P0959R 		PM957R 
 40 mm 1/1	PO188R KIEFER Atraumatic universal and lung forceps	10 mm / 310 mm	P0731R 	PM975R 	P0959R 	PM973R  Ø 5 mm	
 62 mm 1/1	PO215R Intestinal grasping forceps, jaw length 60 mm	10 mm / 310 mm	P0768R 	PM974R 	P0959R 		

Grasping Forceps

AdTec® Monopolar


Grasping Forceps



	Complete Instrument		Single Parts			Individual Components
			Jaw insert	Insulated outer tube	Handle	
 28 mm 1/1	PO186R CROCE/OLMI Combination fixation and dissection forceps, especially for gall bladder	5 mm / 310 mm	 PO729R	 PM973R	 P0959R	 PM975R Ø 5/10 mm  P0958R
 34 mm 1/1	PO105R HEYWOOD-SMITH Tissue- and swab forceps	10 mm / 310 mm	 PO611R	 PM975R	 P0959R	 PM973R Ø 5 mm  P0959R
 41 mm 1/1	PENNINGTON Tissue-grasping forceps Instrument with protruding jaw Parts	10 mm / 370 mm	 PM679R	 PM977R	 P0959R	 P0959R
 19 mm 1/1	PO211R Fixation forceps, medium size serration	5 mm / 310 mm	 PO764R	 PM973R	 P0958R	 PM975R Ø 5/10 mm  PM953R
 19 mm 1/1	PO208R Fixation forceps, medium size serration	5 mm / 420 mm	 PO761R	 PM981R	 P0958R	 PM954R
 18 mm 1/1	PO149R Fixation forceps, jaw Parts fenestrated, fine serration	5 mm / 310 mm	 PO680R	 PM973R	 P0958R	 PM975R Ø 5/10 mm
 18 mm 1/1	PO206R Fixation forceps, jaw Parts fenestrated, medium size serration	5 mm / 310 mm	 PO758R	 PM973R	 P0958R	 PM975R Ø 5/10 mm
 29 mm 1/1	PO130R Atraumatic grasping forceps, jaw part fenestrated	5 mm / 310 mm	 PO646R	 PM973R	 P0958R	 PM975R Ø 5/10 mm  PM957R

AdTec® Monopolar

Grasping Forceps






	Complete Instrument		Single Parts			Individual Components	
			Jaw insert	Insulated outer tube	Handle		
PO129R 18 mm 1/1 Tube fixation and dissecting forceps	5 mm / 310 mm	P0645R PM973R P0958R PM975R Ø 5/10 mm P0958R					
PO200R 18 mm 1/1 Tube fixation and dissecting forceps	5 mm / 420 mm	P0748R PM981R P0958R P0959R					
PO210R 14 mm 1/1 Tube fixation forceps	5 mm / 310 mm	P0763R PM973R P0958R PM975R Ø 5/10 mm P0959R					
PO207R 14 mm 1/1 Tube fixation forceps	5 mm / 420 mm	PM981R P0958R PM953R					
46 mm 1/1 Endo fixation forceps, especially to grasp around the intestine Instrument with protruding jaw Parts	10 mm / 310 mm	PM681R PM974R P0959R PM954R					
PO189R 53 mm 1/1 MASSOUD Delicate forceps for lateral fixation of structures and vascular grafts during dissection and anastomosis For use with flexible 13 mm trocars only	5 mm / 310 mm	P0734R PM973R P0958R PM975R Ø 5/10 mm PM957R					

Grasping Forceps

AdTec® Monopolar


Grasping Forceps



	Complete Instrument		Single Parts			Individual Components	
			Jaw insert	Insulated outer tube	Handle		
 43 mm 1/1	DE BAKEY Atraumatic clamp	5 mm / 310 mm	PM657R	PM973R	P0959R	PM975R Ø 5/10 mm	P0958R
 53 mm 1/1	DE BAKEY Atraumatic clamp	5 mm / 420 mm	P0657R	PM981R	P0959R		P0959R
 68 mm 1/1	DE BAKEY Atraumatic clamp	5 mm / 310 mm	PM658R	PM973R	P0959R	PM975R Ø 5/10 mm	
 35 mm 1/1	DE BAKEY Atraumatic clamp	10 mm / 370 mm	PM651R	PM977R	P0959R		PM953R
 55 mm 1/1	DE BAKEY Atraumatic clamp	10 mm / 220 mm	PM712R	PM992R	PM957R		PM954R
 55 mm 1/1	DE BAKEY Atraumatic clamp	10 mm / 370 mm	P0624R	PM977R	P0959R		
 104 mm 1/1	SIMON Renal pole clamp for local ischaemia during tumor resection	10 mm / 310 mm	P0781R	PM974R	P0959R		PM957R

AdTec® Monopolar


Grasping Forceps








	Complete Instrument		Single Parts			Individual Components	
			Jaw insert	Insulated outer tube	Handle		
 43 mm 1/1	CRAFOORD Endo-vascular clamp	5 mm / 310 mm	PM713R	PM973R	PM957R	PM975R Ø 5/10 mm	P0958R
 59 mm 1/1	CRAFOORD Endo-vascular clamp	5 mm / 310 mm	PM719R	PM973R	PM957R	PM975R Ø 5/10 mm	P0959R
 59 mm 1/1	CRAFOORD Endo-vascular clamp	5 mm / 370 mm	P0684R	PM976R	P0959R	PM978R Ø 5/10 mm	
 38 mm 44 mm 1/1	SATINSKY Atraumatic clamp	5 mm / 310 mm	PM711R	PM973R	PM957R	PM975R Ø 5/10 mm	PM953R
 51 mm 56 mm 1/1	SATINSKY Atraumatic clamp	5 mm / 220 mm	PM690R	PM970R	PM957R	PM971R Ø 5/10 mm	PM954R
 51 mm 56 mm 1/1	SATINSKY Atraumatic clamp	5 mm / 310 mm	PM616R	PM973R	PM957R	PM975R Ø 5/10 mm	PM957R

AdTec® Monopolar

Grasping Forceps





	Complete Instrument		Single Parts			Individual Components	
			Jaw insert	Insulated outer tube	Handle		
 49 mm 1/1	DE BAKEY Pulmonary clamp acc. to LINDER, s-shaped For use with flexible 13 mm trocars only	5 mm / 220 mm	PM693R 	PM970R 	PM957R 	PM971R  Ø 5/10 mm	P0958R 
 49 mm 1/1	DE BAKEY Atraumatic fixation forceps acc. to LINDER, s-shaped For use with flexible 13 mm trocars only	5 mm / 310 mm	PM689R 	PM973R 	PM957R 	PM975R  Ø 5/10 mm	P0959R 
 49 mm 1/1	DE BAKEY Atraumatic fixation forceps acc. to LINDER, s-shaped For use with flexible 13 mm trocars only	5 mm / 310 mm	PM671R 	PM973R 	PM957R 	PM975R  Ø 5/10 mm	PM953R  PM954R 



AdTec® MONOPOLAR

Grasping Forceps













	Complete Instrument		Single Parts			Individual Components	
			Jaw insert	Insulated outer tube	Handle		
 17 mm 1/1	Pike-mouth forceps with reservoir in the jaws	5 mm / 220 mm	P0634R	PM970R	P0959R	PM971R	P0958R
 17 mm 1/1	Pike-mouth forceps with reservoir in the jaws	5 mm / 310 mm	P0613R	PM973R	P0959R	PM975R	P0959R
 18 mm 1/1	Large pike-mouth forceps extra-wide jaws	10 mm / 310 mm	PM621R	PM974R	P0959R		
 33 mm 1/1	Fixation forceps, 2 x 3 teeth	10 mm / 310 mm	P0623R	PM974R	P0959R		PM953R
 14 mm 1/1	Grasping forceps, 2 x 4 teeth	5 mm / 310 mm	P0648R	PM973R	P0959R	PM975R	PM954R
 14 mm 1/1	Grasping forceps, 2 x 4 teeth	5 mm / 420 mm	P0749R	PM981R	P0959R		
 18 mm 1/1	Grasping forceps, 1 x 2 teeth	5 mm / 310 mm	PM649R	PM973R	P0959R	PM975R	PM957R
 15 mm 1/1	Grasping forceps, fine teeth	5 mm / 310 mm	P0650R	PM973R	P0959R	PM975R	


Grasping Forceps

AdTec® Monopolar

Grasping Forceps




Complete Instrument		Single Parts				Individual Components	
			Jaw insert	Insulated outer tube	Handle		
 19 mm 1/1	P0111R NELSON Grasping forceps, sharp	5 mm / 310 mm	 P0617R	 PM973R	 P0959R	 PM975R Ø 5/10 mm	 P0958R
 33 mm 1/1	P0224R Tenaculum Grasping forceps	5 mm / 310 mm	 P0777R	 PM973R	 P0959R	 PM975R Ø 5/10 mm	 P0959R
 42 mm 1/1	P0225R Tenaculum Grasping forceps	10 mm / 310 mm	 P0778R	 PM974R	 P0959R		 P0959R
 23 mm 1/1	P0194R Lymph node grasping forceps	5 mm / 220 mm	 P0739R	 PM970R	 P0959R	 PM971R Ø 5/10 mm	 PM953R
 36 mm 1/1	Endo lymph node fixation forceps	10 mm / 310 mm	 PM688R	 PM974R	 P0959R		 PM954R




Grasping Forceps

AdTec® Monopolar

Biopsy Forceps




	Complete Instrument	Single Parts			Individual Components
		Jaw insert	Insulated outer tube	Handle	
 31 mm 1/1	Large biopsy forceps	10 mm / 310 mm	P0705R	PM974R	P0958R
 15 mm 1/1	Sample excision forceps	5 mm / 310 mm	P0626R	PM973R	P0958R PM975R Ø 5/10 mm
 15 mm 1/1	Sample excision forceps	5 mm / 420 mm	P0750R	PM981R	P0958R
 14 mm 1/1	Biopsy forceps with teeth	5 mm / 310 mm	P0652R	PM973R	P0958R PM975R Ø 5/10 mm
 14 mm 1/1	Biopsy forceps with teeth	5 mm / 420 mm	P0756R	PM981R	P0958R


Biopsy Forceps

AdTec® Monopolar

Biopsy Forceps






	Complete Instrument		Single Parts			Individual Components
			Jaw insert	Insulated outer tube	Handle	
 14 mm 1/1	Biopsy forceps cutting	5 mm / 220 mm	PM635R	PM970R	P0958R	PM971R Ø 5/10 mm
 14 mm 1/1	Biopsy forceps cutting	5 mm / 310 mm	P0627R	PM973R	P0958R	PM975R Ø 5/10 mm
 14 mm 1/1	Biopsy forceps cutting	5 mm / 420 mm	P0759R	PM981R	P0958R	
 14 mm 1/1	FRANGENHEIM Biopsy forceps cutting	5 mm / 310 mm	PM653R	PM973R	P0958R	PM975R Ø 5/10 mm
 14 mm 1/1	FRANGENHEIM Biopsy forceps cutting	5 mm / 420 mm	P0751R	PM981R	P0958R	








AdTec® Monopolar

Needle Holders








	Complete Instrument	Single Parts			Individual Components
		Jaw insert	Insulated outer tube	Handle	
 PO550R 1/1	DUROGRIP Needle holder, with tungsten carbide inserts 5 mm / 310 mm	PO628R 	PM973R 	P0959R 	PM975R  Ø 5/10 mm
 PO559R 1/1	DUROGRIP Needle holder, with tungsten carbide inserts, curved 5 mm / 310 mm	PO668R 	PM973R 	P0959R 	PM975R  Ø 5/10 mm
 1/1	Knot pushing forceps 5 mm / 310 mm	PM654R 	PM973R 	P0958R 	PM975R  Ø 5/10 mm



AdTec® Monopolar


Thorax Instruments








	Complete Instrument		Single Parts			Individual Components
			Jaw insert	Insulated outer tube	Handle	
 1/1	METZENBAUM Scissors, serrated, tips curved to the left, with carbide inserts, distal shaft bent For use with flexible 7 or 13 mm trocars only	5 mm / 310 mm	P0636R + PM940R 	PM983P 	P0958R 	 
 1/1	Peritoneal scissors , serrated, distal shaft bent For use with flexible 7 or 13 mm trocars only	5 mm / 310 mm	PM637R + PM940R 	PM983P 	P0958R 	
 1/1	MARYLAND Fixation and dissecting forceps, curved, distal shaft bent For use with flexible 7 or 13 mm trocars only	5 mm / 310 mm	P0669R + PM940R 	PM983P 	P0958R 	
 1/1	HEYWOOD-SMITH Tissue and swab forceps, distal shaft bent For use with flexible 13 mm trocars only	5 mm / 310 mm	P0641R + PM940R 	PM983P 	P0959R 	

AdTec® Monopolar

Thorax Instruments




Complete Instrument		Single Parts			Individual Components	
		Jaw insert	Insulated outer tube	Handle		
	Pike-mouth forceps with reservoir in the jaws, distal shaft bent For use with flexible 7 or 13 mm trocars only 1/1	5 mm / 310 mm P0642R + PM940R 	PM983P 	P0959R 	P0958R 	
	NELSON Grasping forceps, sharp, distal shaft bent For use with flexible 7 or 13 mm trocars only 1/1	5 mm / 310 mm P0638R + PM940R 	PM983P 	P0959R 	P0959R 	
	Sample excision forceps, distal shaft bent For use with flexible 7 or 13 mm trocars only 1/1	5 mm / 310 mm PM643R + PM940R 	PM983P 	P0959R 	PM953R 	
					PM954R 	
					PM957R 	

Monopolar Laparoscopic Electrodes



The Modular Principle Also Includes Important Safety Aspects

- modular design simplifies disassembly, cleaning, maintenance and sterile preparation
- bayonet connection
- Individual combinations possible
- fast exchange of Individual Components due to modular design
- ceramic tips
- insulation up to the working tip



GK383R		Hook electrode, J-shaped
GK384R		Hook electrode, L-shaped
GK385R		Ball electrode
GK386R		Spatula electrode
GK393R		Hook electrode, J-shaped
GK394R		Hook electrode, J-shaped, 90° curved to the right
GK395R		Ball electrode, curved

MONOPOLAR LAPAROSCOPIC ELECTRODES


GK390R	
GK389R	Pin electrode only (pack of 5 pcs.) product is to be reprocessed only once (delivery unsterile)
GK373R	 Handle with inner tube, diam. 5 mm
GK371P	 Separate outer tube, diam. 5 mm









Art. no.	Monopolar connecting cable 3.5 m; on the Instrument side for a 4 mm contact pin
GN202	AESCLAP®, ACMI, Berchtold, Martin
GK239	ERBE ICC / ACC, KARL STORZ
GK245	ERBE T
GK246	Valleylab, ACMI, Conmed, Eschmann, Erbe ICC Intern., US-Norm
GK299	 the electrodes

AdTec® Monopolar

Spare Parts

		
P0950400	AESCLAP®, Wolf, Storz (Click-line), Olympus (HiQ+), including cap, package of 5 pieces	 P0950R
P0950410	Storz, old type, Package of 5 pieces	
		
P0950420	Olympus (HiQ) diam. 3 mm, package of 5 pieces	 P0951R

Art. no.	HF contact pin for AdTec® Monopolar handle PM950R, PM951R, PM953R, PM955R only	
		
PL170251	AESCLAP®, Wolf, Storz (Click-line), Olumpus (HiQ+), including cap, package of 5 pieces	 PM950R
PL170256	Olympus (HiQ), diam. 3 mm, package of 5 pieces	 PM951R
PL170258	Storz, old type, package of 5 pieces	 PM953R
		
PL170250	Luer-Lock adapter replaces HF-pin for rinsing the Instrument, package of 5 pieces for handles P0950R, P0951R, PM950R, PM951R, PM953R, PM955R	
		 PM955R

		
PM950250	Rotation knob, package of 5 pieces	 PM950R
		 PM951R

PM630201	220 mm
PM600201	310 mm
PM618201	370 mm
PM434201	420 mm




Art. no.	Inner tube for thorax Instruments (PM940R), consisting of:
----------	--


	
PM009250	Inner tube, curved on flat, for Instruments, length 310 mm

	
--	---


PM009251	Flexible tension rod	
PM970R	220 mm	5/5 mm
PM971R	220 mm	5/10 mm
PM992R	220 mm	10/10 mm
PM973R	310 mm	5/5 mm
PM975R	310 mm	5/10 mm
PM974R	310 mm	10/10 mm
PM976R	370 mm	5/5 mm
PM978R	370 mm	5/10 mm
PM977R	370 mm	10/10 mm
PM981R	420 mm	5/5 mm
PM982R	420 mm	10/10 mm
PM983P	310 mm	5/5 mm




Endo Rack



- Universal silicone storage rack for Instruments
- Variable Instrument storage
- Fold down handles for easier Instrument removal
- Optional trocar storage




Art. no.	Storage rack set for MIS Instruments, standard
	
	consisting of: Storage rack, standard, perforated tray and silicon mat for Instruments with shaft length up to 370 mm Dimensions: 540 x 253 x 166 mm
PL960R	

For more information see brochure number C22111.



JK444 + JK489

JK444 bottom and JK489 lid
 External dimensions: 592 x 285 x 205 mm

Art. no.	Storage rack set for MIS Instruments, long
	
	consisting of: Storage rack, long, perforated tray and silicon mat for Instruments with shaft length up to 420 mm Dimensions: 660 x 253 x 166 mm
PL961R	



JN445 + JK490


JN445 bottom and JK490 lid
 External dimensions: 711 x 285 x 205 mm
 NB: The Endocontainer can only be sterilized in sterilizers with the appropriate chamber size (at least 720 mm working depth).

AdTec® Monopolar


Standard Instruments

Disassembly


Assembly




I Deactivate ratchet on ratchet handle
Flip down red toggle switch. **N.B.: for P0959R only**




II Disconnect shaft and handle
Pull rotation knob back until it locks in place. A red ring becomes visible.
Important: Do not hold on to the movable handle component.




III Assemble shaft I
Slide the metal inner tube onto the jaw insert with its four spring tongues towards the working end until the tongues snap in place.




IV Assemble shaft II
Slide the outer insulation tube as far as it will go over the inner metal tube.



V Remove shaft from handle
Remove the shaft by pulling the Instrument shaft out of the handle.




VI Disassemble shaft I
Hold the jaw insert at the tip and remove the outer insulation tube.




VII Prepare handle
If the red ring is not visible, pull the rotation knob back. Deactivate the ratchet by flipping down the red toggle switch. **N.B.: for P0959R only**



VIII Connect shaft and handle
Hold the shaft firmly by its tip and the handle behind the rotation knob. Push the shaft into the handle until it locks in place and the red ring disappears.
Important: Do not hold on to the movable handle component.



IX Disassemble shaft II
Hold the jaw insert at the tip and remove the inner metal tube.



X Function test
Test the assembly of the Instrument by opening and closing the handle to ensure the jaw is opening and closing.

AdTec® Monopolar

Curved Thoracic Instrument



Disassembly


Assembly




I Pull back the spring mechanism indicated with an arrow as far as it will go. Hold at this position and at the same time pull out the shaft. **Note:** In the case of handle with a lock (PM954R, PM957R), open the ratchet Completely.




II Pull off the insulated, flexible outer tube from the shaft.




III Push the flexible tension rod into the inner tube.




IV Push the tension rod forwards until the receptacle of the jaw part (screw thread) is clearly visible.




V Push flexible tension rod towards the jaw part (approx. 1 cm).




VI Unscrew the jaw part as indicated by the arrow.




VII Screw on jaw part. Set the position of the jaw part relative to the shaft to meet the surgeon's needs. To do this, the thread connection must be slightly slackened (max. 1 turn).



VIII Push jaw part with the flexible tension rod into the inner tube with the 4 slots.




IX After unscrewing the jaw part, pull the tension rod backwards.



X Remove the Complete tension rod in a backward direction from the inner tube.




XI Push insulated flexible outer tube PM983P onto the inner tube.



XII Press movable ring part of the handle forwards and hold at this position. Pull back spring mechanism of the axial handle as far as it will go and hold it. Take the shaft at the distal tip, push shaft into the spring mechanism as far as it will go and join with a slightly rotary movement.

AdTec® Monopolar

Care



Lubricate the defined places on the Instruments (Instruments jaw and red toggle switch) with Aesculap® Sterilit Oil after each cleaning cycle.



Art. no.	AESCULAP® Oil Spray
JG600	AESCULAP® Oil Spray in aerosol can, for maintenance of Instruments prior to sterilization.
JG598	AESCULAP® Instrument Oil for box locks, ratchets etc. prior to sterilization.



GK469R	Cleaning brush diam. 3.7 mm, length 550 mm for inner tube and outer tube diam. 5 mm. Package of 6 pieces.
--------	---



PM995R	Cleaning brush diam. 10 mm, length 550 mm for outer tube diam. 5 mm and diam. 10 mm. Package of 6 pieces.
--------	---

AdTec® Monopolar

Numerical Index

Catalogue Number	Page	Catalogue Number	Page	Catalogue Number	Page	Catalogue Number	Page	Catalogue Number	Page	Catalogue Number	Page
GK239	57	JK490	61	PM630201	59	PM701R	17	P0001R	10	PO114R	26
GK245	57	JN445	61	PM630R	11	PM710R	9	P0002R	10	PO117R	40
GK246	57			PM635R	47	PM711R	37	P0003R	10	PO118R	34
GK299	57	PL170250	58	PM637R	51	PM712R	35	P0004R	6	PO127R	20
GK370P	55	PL170251	58	PM640R	27	PM713R	37	P0006R	10	PO128R	22
GK371P	56	PL170256	58	PM643R	53	PM719R	37	P0011R	8	PO129R	32
GK372R	55	PL170258	58	PM647R	15	PM940R	51, 53	P0014R	10	PO130R	30
GK373R	55 f.	PL960R	61	PM649R	41	PM950250	58	P0022R	6	PO132R	40
GK374R	55	PL961R	61	PM651R	35	PM950R	58	P0023R	8	PO134R	40
GK375R	55	PM009250	59	PM653R	47	PM951R	58	P0025R	10	PO137R	34
GK376P	55	PM009251	59	PM654R	49	PM953R	7 ff.	P0027R	6	PO141R	16
GK383R	55	PM171R	36	PM657R	35	PM954R	7 ff.	P0028R	8	PO147R	18
GK384R	55	PM173R	36	PM658R	35	PM955R	58	P0030R	8	PO149R	30
GK385R	55	PM175R	38	PM670R	29	PM957R	7 ff.	P0031R	8	PO151R	18
GK386R	55	PM176R	38	PM671R	39	PM970R	7 ff.	P0090R	6	PO152R	36
GK389R	56	PM177R	36	PM679R	31	PM971R	7 ff.	P0091R	6	PO185R	24
GK390R	56	PM178R	36	PM681R	33	PM973R	7 ff.	PO100R	14	PO186R	30
GK393R	55	PM434201	59	PM687R	29	PM974R	7 ff.	PO101R	18	PO187R	14
GK394R	55	PM600201	59	PM688R	43	PM975R	7 ff.	PO102R	12	PO188R	28
GK395R	55	PM604R	7	PM689R	39	PM976R	17 ff.	PO103R	16	PO189R	32
GK469R	67	PM609R	29	PM690R	37	PM977R	27 ff.	PO104R	28	PO190R	18
GN202	57	PM612R	29	PM693R	39	PM978R	17 ff.	PO105R	30	PO191R	24
		PM614R	15	PM696R	29	PM981R	7 ff.	PO107R	40	PO194R	42
JG598	67	PM616R	37	PM697R	7	PM982R	27, 59	PO109R	20	PO195R	14
JG600	67	PM618201	59	PM698R	9	PM983P	51 ff.	PO111R	42	PO196R	26
JK444	61	PM618R	25	PM699R	13	PM992R	11, 35, 59	PO112R	24	PO197R	26
JK489	61	PM621R	41	PM700R	17	PM995R	67	PO113R	26	PO199R	18

AdTec® Monopolar

Numerical Index

Catalogue Number	Page	Catalogue Number	Page	Catalogue Number	Page	Catalogue Number	Page	Catalogue Number	Page	Catalogue Number	Page
P0200R	32	P0229R	22	P0608R	13	P0648R	41	P0739R	43	P0769R	23
P0201R	40	P0230R	20	P0610R	29	P0650R	41	P0740R	15	P0770R	23
P0202R	12	P0231R	14	P0611R	31	P0652R	45	P0741R	27	P0771R	17
P0203R	18	P0232R	24	P0613R	41	P0655R	9	P0742R	27	P0772R	17
P0204R	24	P0233R	22	P0615R	21	P0657R	35	P0744R	11	P0773R	17
P0205R	26	P0234R	22	P0617R	43	P0661R	17	P0746R	7	P0774R	21
P0206R	30	P0235R	22	P0618R	25	P0667R	19	P0747R	19	P0775R	25
P0207R	32	P0236R	22	P0619R	27	P0668R	49	P0748R	15, 33	P0776R	25
P0208R	30	P0237R	12	P0620R	27	P0669R	51	P0749R	41	P0777R	43
P0209R	26	P0238R	12	P0623R	41	P0678R	19	P0750R	45	P0778R	43
P0210R	32	P0500R	44	P0624R	35	P0680R	31	P0751R	47	P0781R	35
P0211R	30	P0501R	44	P0625R	11	P0682R	11	P0752R	13	P0782R	23
P0212R	26	P0502R	46	P0626R	45	P0684R	37	P0753R	19	P0783R	21
P0213R	16	P0505R	44	P0627R	47	P0705R	45	P0754R	9	P0785R	25
P0214R	18	P0509R	44	P0628R	49	P0715R	9	P0755R	25	P0786R	23
P0215R	28	P0510R	46	P0631R	7	P0716R	9	P0756R	45	P0787R	23
P0216R	22	P0511R	44	P0632R	13	P0717R	13	P0757R	27	P0788R	23
P0217R	22	P0512R	46	P0633R	19	P0718R	13	P0758R	31	P0789R	23
P0218R	16	P0550R	48	P0634R	41	P0727R	7	P0759R	47	P0790R	7
P0219R	16	P0559R	48	P0636R	51	P0728R	25	P0761R	31	P0791R	7
P0220R	16	P0600R	11	P0638R	53	P0729R	31	P0762R	27	P0950400	58
P0221R	20	P0601R	11	P0641R	51	P0730R	15	P0763R	33	P0950410	58
P0222R	24	P0602R	11	P0642R	53	P0731R	29	P0764R	31	P0950420	58
P0223R	24	P0603R	7	P0644R	21	P0734R	33	P0765R	27	P0950R	58
P0224R	42	P0605R	15	P0645R	33	P0735R	19	P0766R	17	P0951R	58
P0225R	42	P0606R	17	P0646R	31	P0736R	25	P0767R	19	P0958R	7 ff.
P0228R	34	P0607R	19	P0647R	23	P0738R	9	P0768R	29	P0959R	7 ff.

Notes

Notes

AESCULAP® – a B. Braun brand

B. Braun Australia Pty Ltd | Level 2, 4 Burbank Place, Norwest NSW 2153, Australia |
Tel. 1800 251 705 | info.au@bbraun.com | www.bbraun.com.au
B. Braun New Zealand | Level 8, 139 Quay Street, Auckland, 1010, New Zealand |
Customer Care 0800 227 286 | Fax (09) 373 5601 | www.bbraun.co.nz